



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Rubber Cones For Filtration Flask

Product Image



Description

- Used to avoid glass to glass contact.
- Used on many Filtering Apparatus, Flasks, Buchner funnel etc.
- Available in Natural & Silicon Rubber.
- Available in Pack of 25.

Silicon Rubber Cones:

Natural Rubber Cones:

**Catalog
No.** **Sizes**

**Catalog
No.** **Sizes**

101354-A	250/500ml	101354-E	250/500ml
101354-B	1/2 Litre	101354-F	1/2 Litre
101354-C	5 Litre	101354-G	5 Litre
101354-D	10 Litre	101354-H	10 Litre

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.