

#### Address: JAMBU PERSHAD & SONS

6275/22 Nicholson Road, Ambala Cantt, Haryana, INDIA

Pin: 133001

**Email:** sales@japson.com japsonambala@yahoo.com

Website: www.japson.com Phone:

+91-171-4006897

# **Single Rotary Tablet Press**

## **Product Image**



## **Description**

Our single rotary tablet press is widely used in the pharmaceutical and related industries since it is simple to operate and requires little maintenance. The quick change-over and volume adjustment feature aid in the production of tablets of various sizes and forms.

The machine's body is made of SS 304, while the contact portions are made of SS16, providing it considerable strength and assuring durability. The machine has pre-pressing

and main-pressing mechanisms. The compression force is adjustable up to 6 tons.

High Quality product for use in Pharmaceutical Industry and Laboratory. Following is made in compliance with International Quality Standards and are manufactured using high-grade components and materials.

\_

Single Rotary Tablet Press, 12 Stations	<b>Catalog No.</b> 101277 (A)
Single Rotary Tablet Press, 16 Stations	<b>Catalog No.</b> 101277 (B)
Single Rotary Tablet Press, 20 Stations	<b>Catalog No.</b> 101277 (C)
Single Rotary Tablet Press, 23 Stations	<b>Catalog No.</b> 101277 (D)

#### **Disclaimer**

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.