



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Small Sample Adaptor

Product Image



Description

The Small Sample Adaptor, consisting of a coupling nut,
an extension link, a spindle set
(spindle#21,#27,#28 and #29),
a sample chamber and an adjustable base

Salient Features

- ✦ Small Sample Adaptor is recommended for small sample volume.
- ✦ Small Sample Adaptor widely used for Digital Viscometer
- ✦ 10 to 20 ml sample is needed for test in Small Sample Adaptor.
- ✦ Sample temperature can be controlled by using Refrigerated Circulating Bath (Optional)
- ✦ The Sample Chamber fit into a water jacket.

* Small Sample Adaptor is an Optional Accessory with Digital Viscometer

Note : Small Sample Adaptor order along with Viscometer only. It **cannot** be attached with the main instrument later on.

Small Sample Adaptor

Catalog No. 101127

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.