



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Smoke Cell

## Product Image



## Description

---

For demonstrating Brownian motion. Designed for use with a microscope. A short cylindrical cell, closed by a cover glass, is held vertically in an opaque channel 90x90x30mm. Illumination is by an integral filament bulb & cylindrical glass lens. The cell enables Brownian motion to be observed for at least 15 minutes after charging the cell with smoke. Fitted with 4mm socket terminals for connection with 12V supply.

Catalog No.	Particulars
JZ02960	Smoke Cell

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy. Other configurations and superior quality available for specific use or for developed countries.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.