



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Specimen Polishing Machine

## Product Image

---



## Description

---

Specimen Polishing Machine can be used for coarse & fine grinding, rough polishing & fine polishing. It is the perfect metallographic specimen machine for samples preparation.

Catalog No.	Plate Disc	Diameter of plates	RPM	Remarks	Coolant Facility	Dimensions/ Weight in Inches
7001 PMD	2	200mm	300-650	Input Voltage AC 220Volts / 50Hz 50 Hz single phase Double Motor of 0.5 HP each	Inlet and Outlet Coolant Pipe	48x20x26 / 75Kg
7001 PMS	1	200mm	300-650	Input Voltage AC 220Volts / 50Hz Single Phase Single motor of 0.5 HP	Inlet and Outlet ,Coolant Pipe	29x21x20 / 45Kg

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.