



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Specimen Polishing Machine

Product Image



Description

Specimen Polishing Machine can be used for coarse & fine grinding, rough polishing & fine polishing. It is the perfect metallographic specimen machine for samples preparation.

| Catalog No. | Plate Disc | Diameter of plates | RPM | Remarks | Coolant Facility | Dimensions/ Weight in Inches |
|-------------|------------|--------------------|---------|--|--------------------------------|------------------------------|
| 7001 PMD | 2 | 200mm | 300-650 | Input Voltage AC 220Volts / 50Hz 50 Hz single phase Double Motor of 0.5 HP each | Inlet and Outlet Coolant Pipe | 48x20x26 / 75Kg |
| 7001 PMS | 1 | 200mm | 300-650 | Input Voltage AC 220Volts / 50Hz Single Phase Single motor of 0.5 HP | Inlet and Outlet ,Coolant Pipe | 29x21x20 / 45Kg |

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.