



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Stability Chamber GMP MODEL

Product Image



Description

This equipment provides a long-time stable environment of temperature, humidity and illumination for the Pharmaceutical industry medicine invalidation evaluation. It's applied in the medicine accelerating test, growth test, high humidity test and high illuminance test.

Features:

The Inner chamber and Outer chamber; 304 high-quality stainless steel, full arc design, easy to clean and maintain.

The Unique Air Duct design ensures the uniform temperature and humidity in the box.

Remote monitoring via mobile phones and computers.

It has three levels of Password Authority Management function.

Data recording system can record more than 1 million pieces of data, even if the power failure, the data will not be lost, can be exported through the U disk.

Conforms to the national pharmacopoeia, GMP, FDA, ICH and other relevant standards for long-term test, intermediate test, accelerated test, with stable and reliable performance to meet the requirements of temperature, humidity, light and other environmental conditions in multiple drug tests.

Cat No.	Description	Internal Size (WxDxH) mm
101504A	STABILITY CHAMBER, 250L	650 X 465 X 830
101504B	STABILITY CHAMBER, 450L	660 x 655 x 1050
101504C	STABILITY CHAMBER, 1000L	800 x 800 x 1580

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.