



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Stefan Constant Apparatus Experiment Apparatus

Product Image



Description

Objective : Stefan Constant apparatus has been designed to study the Stefans law of radiation.

Features : Continuously variable, over load & short circuit protected DC regulated power supply of 0 to 20V. Voltmeter & Current meter are mounted on the front panel to measure the voltage & current. Bulb holder with bulbs 19 Volt is mounted on the front panel.

Catalog No. Particulars

JE21620	Stefan Constant Apparatus Experiment Apparatus
---------	--

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.