



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Stefan Constant Apparatus Experiment Apparatus

Product Image



Description

Objective : Stefan Constant apparatus has been designed to study the Stefans law of radiation.

Features : Continuously variable, over load & short circuit protected DC regulated power supply of 0 to 20V. Voltmeter & Current meter are mounted on the front panel to measure the voltage & current. Bulb holder with bulbs 19 Volt is mounted on the front panel.

Catalog No. Particulars

JE21620 Stefan Constant Apparatus Experiment Apparatus

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.