



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Submarine Gel Electrophoresis, Analyst Mechanism

Product Image



Description

This is a laboratory scale mechanism. In this mechanism we can run 3 different gels & screen up to 32 samples at a time. This elementary unit consists of detachable platinum electrode assembly(2 nos.), perforated lid for vapors, gel trays 13x13cm., 13x6cm.(2 nos.), 6x6.5cm(4nos.), 15well comb 1.5mm, 8-well comb 1.5mm(2 nos.), 2 well comb 3mm, comb stand connecting cord & instruction manual.

Particulars

Catalog No.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.