



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Submarine Gel Electrophoresis, Master Mechanism

Product Image



Description

Submarine master electrophoresis is designed for running long gels. Medium gels can also be accomplished very easily in this mechanism. The elementary unit consists of detachable platinum electrode assembly (2 nos.), perforated lid for vapours, gel trays 25x 13cm., 15 well comb 1.5mm, 3 well comb 3.0mm, 3 well comb 6.0mm, comb stand, connecting cord & instruction manual.

Particulars

Submarine Gel Electrophoresis, Master
Mechanism

Catalog No.

JA20750

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.