



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Test Tube Stand, Plastic

Product Image



Description

Material: Polypropylene

The autovlable stand is perhaps the most commonly used Test tube stand in various laboratories .For better visibility the six place rack is end supported with two columns.The top plates has two holes for holding 25 mm dia test ube,4 holes to hold 16 mm dia tubes whereas the base has corrsponding hemisephericall wells.It also has 6 vertical pins for drying tubes

Catalog No.	Particulars (mm)
JW08401	16 & 25
JW08402	25

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.