



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# Thermometer, Laboratory, Superior

## Product Image



## Description

---

Glass capillary with bulb type, yellow back with printed scale. Length 300mm with blue plastic ring on top.

Different types available :

Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top

Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top

### Catalog No.Particulars

JD14680	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14681	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14682	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14683	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14684	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -50 to + 50 C, 1.0
JD14685	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -50 to + 100 C, 1.0
JD14686	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -100 to + 50 C, 1.0
JD14687	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, -100 to + 100 C, 1.0
JD14688	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, 0 to 220 F, 1.0
JD14689	Thermometer, Laboratory, Superior, Alcohol in plastic case with blue plastic ring on top, 0 to 400 F, 2.0
JD14690	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 50 C, 0.5
JD14691	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 110 C, 1.0
JD14692	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 110 C, 0.5
JD14693	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 150 C, 1.0
JD14694	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 200 C, 1.0
JD14695	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 250 C, 1.0

JD14696	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 300 C, 1.0
JD14697	Thermometer, Laboratory, Superior, Mercury in plastic case with blue plastic ring on top, -10 to 360 C, 2.0

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.