



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Thin Layer Chromatography Apparatus, Brass

Product Image



Description

To Introduce the principles of absorption chromatography using the thin layer technique. This mechanism is having moveable applicator with inbuilt thickness arrangement between 0 to 2 mm (Minimum Div. 0.25 mm). It further consists of the following components. a) Spreader (Applicator) made of electroplated BRASS. b) Perspex base of 114 x 23 cm. to support glass plates. c) Plate rack aluminium, anodised, for ten 20 x

20 cm. or 20 x 10 cm. plates. d) Developing tank with lid. e) Spotting template made of perspex. f) T.L.C. Plates set of five 20 x 20 cm. and two 20 x 5 cm. or set of ten 20 x 10 cm. and two 20 x 5 cm. g) Glass sprayer with rubber bellow, Cap. 100 ml. h) Micro-pipette. i) Scriber for marking lines made of stainless steel. j) Instruction Manual.

Particulars	Catalog No.
Thin Layer Chromatography Apparatus, BrassThin Layer Chromatography Apparatus, Brass	JA20650

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.