



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tissue Culture Dish

Product Image



Description

Tissue Culture Dishes are ideal containers for cell growth and yield. These dishes are available either non-treated or surface treated to satisfy different applications. The uniform surface treatment makes the hydrophobic surface hydrophilic which improves the adhesion of cells and binding of proteins to surface.

Type A Surface Treated, Sterilized

Type B Non-Treated, Sterilized

Type A: Surface Treated, Sterilized

Catalog No.	Diameter (mm)
JW07861	35
JW07862	70
JW07863	90
JW07864	150

Type B: Non-Treated, Sterilized

Catalog No.	Diameter (mm)
JW07865	35
JW07866	70
JW07867	90
JW07868	150

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.