



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tissue Culture Flask

Product Image



Description

Tissue Culture Flasks are extremely perfect for cell growth and yeild. These flasks are available either non-treated or surface treated to satisfy different application. The uniform surface treatment makes the hydrophobic surface hydrophilic which improves the adhesion of cells and binding of proteins to surface. These flasks are available with filter cap or seal cap.

Type A Surface Treated, Sterilized

Type B Non-Treated, Sterilized

Type A: Surface Treated, Sterilized

Catalog No.	Capacity (ml)
JW07841	50
JW07842	250
JW07843	600

Type B: Non-Treated, Sterilized

Catalog No.	Capacity (ml)
JW07844	50
JW07845	250
JW07846	600

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.