



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Tissue Culture Flask with Filter Cap

## Product Image

---



## Description

---

Tissue Culture Flask with Filter Cap has a hydrophobic filter which allows optimal gas exchange and reduce chance of contamination.

**Type A** Surface Treated, Sterilized

**Type B** Non-Treated, Sterilized

### **Type A: Surface Treated, Sterilized**

<b>Catalog No.</b>	<b>Capacity (ml)</b>
JW07851	50
JW07852	250
JW07853	600

### **Type B: Non-Treated, Sterilized**

<b>Catalog No.</b>	<b>Capacity (ml)</b>
JW07854	50
JW07855	250
JW07856	600

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.