



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Tissue Culture Rack

Product Image



Description

Construction: Tissue racks are made of tubular mild steel pipes duly painted for longer life/ made of aluminium.

Racks: Each rack of size 48"x12" (depth) is covered with thick unbreakable plexi glass/sunmica paste duro board, which are terminate and insect proof fire and weather resistant. Distance between each rack is 15".

Lighting: Four fluorescent tubes are properly illuminated by diffused, glare free

fluorescent lights are built in upper portion in worktable. UV Germicidal Light are fitted inside the rack/ work area. All fluorescent tubes chocks and ultra violet germicidal tube chocks are fitted outside of rack. If desired to avoid heat dissipation inside the rack.

Timer: 0-24 Hours Digital Timers can be provided in tissue culture rack.

Control Panels: All control panels are fitted tops of tissue culture racks including On/Off switch of fluorescent tubes and UV tubes. All control panels are fitted separately in each rack with indicator light.

Easy Mobility: The whole rack is provided with four castor wheels at the bottom for easy mobility to any place.

Unit of four shelves with 16 Fluorescent tubes 4 U.V. lights.

Unit of five shelves with 20 Fluorescent tubes 5 U.V. lights.

Unit of six shelves with 24 Fluorescent tubes 6 U.V. lights.

Tissue Culture Rack

Catalog No. 100398

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.