



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# VDRL Shaker

## Product Image



## Description

---

It is very useful for V.D.R.L. grouping & typing as well as other agglutination tests. It is also used for mixing of solution in small flask, beakers, bottles or similar containers. It has a timer incorporated in the electrical circuit which is being cut off automatically as soon as the setting is over. Platform 30x30 cm. Also available with digital timer.

Particulars	Catalog No.
V.D.R.L. Shaker	JA 23830
V.D.R.L. Shaker with Digital RPM & Timer	JA 23831

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.