



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Vacuum Filtration Apparatus

Product Image



Description

Vacuum filtration apparatus are useful to purify large volume samples efficiently in a easy and fast manner. The filters are ultrasonically weld in the cup of apparatus and bottle has a wide mouth for efficient handling. Both Upper and lower cups have graduations.

Membrane Material: PVDF

Catalog No.	Pore Size (um)	Capacity (ml)
JW13411	0.1	150
JW13412	0.1	250
JW13413	0.1	500
JW13414	0.22	150
JW13415	0.22	250
JW13416	0.22	500
JW13417	0.45	150
JW13418	0.45	250
JW13419	0.45	0.45
JW13420	0.22	150
JW13421	0.22	250
JW13422	0.22	500
JW13423	0.45	150
JW13424	0.45	250
JW13425	0.45	500
JW13426	0.45	150
JW13427	0.45	250
JW13428	0.45	500

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.