



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Vacuum Tray Dryer

Product Image



Description

This Vacuum Tray Dryer enables total moisture vaporization and product drying. It is appropriate for drying thermally sensitive materials that are quickly resolved, polymerized, or degraded at higher temperatures.

The materials to be dried are placed on trays within the vacuum tray dryer, and pressure is decreased using a vacuum pump. Sterilization can be performed before drying, during which no matter is permitted to enter the product. The dryer door is firmly closed, and steam is fed through the gap between the trays and the jacket, allowing heat transmission by conduction. The condenser receives water vapors from the feed and once dried, The vacuum pump is deactivated, and the dried product is collected from the trays.

High Quality product for use in Pharmaceutical Industry and Laboratory. Following is made in compliance with International Quality Standards and are manufactured using high-grade components and materials.

-

Vacuum Tray Dryer, 6 Trays	Catalog No. 101254 (A)
Vacuum Tray Dryer, 12 Trays	Catalog No. 101254 (B)
Vacuum Tray Dryer, 18 Trays	Catalog No. 101254 (C)
Vacuum Tray Dryer, 24 Trays	Catalog No. 101254 (D)
Vacuum Tray Dryer, 36 Trays	Catalog No. 101254 (E)
Vacuum Tray Dryer, 48 Trays	Catalog No. 101254 (F)
Vacuum Tray Dryer, 96 Trays	Catalog No. 101254 (G)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.