



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Viscosity Apparatus Stokes Method

## Product Image

---



## Description

---

A Borosilicate glass tube of about 30 cm diameter and one meter long is made to stand vertically firmly in between two upright bars with the help of clamps on a heavy cast iron base. Both sides of the tube are open. On the lower end a rubber cork without holes is fitted whereas on upper end a rubber cork with a funnel is provided. With two rubber band index mark and 6 steel balls of different diameter. Without lead shot.

Viscosity Apparatus Stokes Method

**Catalog No. 100411**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.