



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Viscosity Apparatus Stokes Method

Product Image



Description

A Borosilicate glass tube of about 30 cm diameter and one meter long is made to stand vertically firmly in between two upright bars with the help of clamps on a heavy cast iron base. Both sides of the tube are open. On the lower end a rubber cork without holes is fitted whereas on upper end a rubber cork with a funnel is provided. With two rubber band index mark and 6 steel balls of different diameter. Without lead shot.

Viscosity Apparatus Stokes Method

Catalog No. 100411

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.