



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Vortex Mixer

Product Image



Description

Salient Features

- ✦ Touch and continuous operations
- ✦ 2500RPM, Fixed Speed Vortex Mixer
- ✦ Used for various mixing applications.
- ✦ Specially designed vacuum suction feet hold to the work
- ✦ Supplied with Standard Head

The Vortex mixer widely used in life science, physical and chemical analysis fields, Analytical & Microbiology laboratories to suspend cell & Mixing reagent to mix sample and dilutant.

SPECIFICATIONS	Model	LMVM20
	Speed Range	2500rpm
	Shaking movement	Orbital
	Orbital diameter	4 mm
	Motor Rating input	58W
	Motor Rating Output	10W
	Run Type	Continuos / Touch Operation
	Ambient Temperature range	5°C to 40°C
	Display	None
	Power	60W
	Power supply	AC 220V / 50Hz
	Dimension LxWxH (mm)	170x120x140
	Weight	2.5kg

Vortex Mixer, Euro

Catalog No. 101123

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.