



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Vortex Shaker (Test Tube Shaker)

Product Image



Description

For shaking liquids in Test Tubes. It is a vibration type with variable speed. It gives short rapid eccentric strokes. Neoprene cap imparts vigorous agitation to tube or small flasks. Shaking starts when vessel is pressed against the cup. With ON/OFF switch, speed control and indicator. Vacuum rubber fitted at the bottom to prevent movement. Also available with adapters for multiple tube.

Particulars	Catalog No.
Vortex Shaker (Test Tube Shaker)	JA 23810
Vortex Shaker (Test Tube Shaker) with Adapter for Multiple Tubes	JA 23811

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.