



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type)

Product Image



Description

Available in various chamber sizes with cover. Three delivery sieves 20, 40 and 60 mesh respectively, each with plastic receiver. Works on 220 volts AC. A removeable glass plate forms the face of the chamber and permits observation during operation.

Particulars

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type) With
Chamber Size 40 mm x 25 mm

Catalog No

JA23850

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type) With Chamber Size 65 mm x 25 mm JA23851

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type) With Chamber Size 40 mm x 25 mm (All critical parts of Stainless Steel) JA23852

Willy Grinder Or Lab Willy Mill (Arthur Thomas Type) With Chamber Size 60 mm x 25 mm (All critical parts of Stainless Steel) JA23853

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.