

Address: JAMBU PERSHAD & SONS

6275/22 Nicholson Road, Ambala Cantt, Haryana, INDIA

Pin: 133001

Email:

sales@japson.com japsonambala@yahoo.com Website:

www.japson.com

Phone:

+91-171-4006897

Zener Diode Characteristics Apparatus

Product Image



Description

Zener diode characteristics apparatus (forward & reverse v-i characteristics)

Objective: to study forward & reverse v-i characteristics of zener diode.

Features: instrument comprises of dc regulated power supply 0-15vdc/150ma, two round meters for voltage & current measurement, 3 zener diodes mounted behind the panel, connections of supplies, meters & zener diode brought out at 4mm sockets.

Catalog No. Particulars

JE20020 Zener Diode Characteristics Apparatus

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.