



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Anaerobic Culture Jar 3.5 Litre with Vacuum cum Pressure Gauge

Product Image



Description

Transparent anaerobic culture jar of 3.5 litre capacity which is made of polycarbonate. The lid of the jar is made of aluminium pressure die casting fitted with pressure valve, vacuum valve, safety valve & compound pressure cum vacuum gauge, three finger aluminium clamp and sealing "O" ring and air/gas chuck.

This is recommended only where inside condition of air of jar is removed thru vacuum pump and CO₂ is filled in the jar thru CO₂ gas cylinder to generate anaerobic conducting in the jar. Also supplied with plain aluminium pressure die casted lid which can be fitted when readymade gas packs are used to generate the anaerobic conditions inside the jar.

Supply completed with following accessories :

- Petri plate carrier, stainless steel, for 10 petri dishes.
- Test tube carrier, stainless steel.

Optional consumables at extra cost:

- Anaerobic Indicator, Box of 10 Tablets.
- Anaerobic Gas Pack, Price per Pkt. of 5 Sachets.

Spare Parts at extra cost:

- Silicon ring
- Petri Plate Carrier, Stainless steel, for 10 Petri Dishes.
- Test Tube Carrier, Stainless Steel.
-

Particulars

Anaerobic Culture Jar 3.5 Litre with Vacuum cum Pressure Gauge

Catalog No.

JA90350

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.