



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Analog Signal Sampling and Reconstruction Kit Experiment Apparatus

Product Image



Description

AF Signal generator 20 Hz to 20 KHz, 5 V p-p
Clock generator 2 KHz to 32 KHz to study:
Sample and hold
Natural sampling
Flat top sampling
Sampling clock width control 10% to 90
Low pass filter
Power supply +/- 12 V, 500 mA and +/- 5 V, 50 mA
Required numbers of patch cords and operating manual

Catalog No.	Particulars
--------------------	--------------------

JE21840	Analog Signal Sampling and Reconstruction Kit Experiment Apparatus
---------	---

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.