



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Aneroid Sphygmomanometer

Product Image



Description

Round shaped dial diameter 147mm, special precise mechanism, metal castors, 360 degree swivel for easy reading, gradient adjustable for easy reading, 4-wheel adjustable height, metal basket, nylon cuff with D-ring, 2-tube adult latex bladder, standard latex bulb, air-release valve with spring, standard end valve, spiral PVC tube (1.8 m) and plastic connector.

Catalog No.	Particulars
JH10550	Aneroid Sphygmomanometer

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.