



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Animal House

Product Image



Description

RABBIT CAGE :

With Heavy gauge mesh mounted on a sturdy frame of Mild Steel, angle size 18"x14"x14" with door in front. 2" Deep Tray, Draw of G.I. Sheet is provided with slides in an angle attached to the sides of the cage and is removable for cleaning. Finished with Aluminum Lacquer. Card holder is provided along.

RABBIT CAGE Stainless Steel Accessories :

Complete with all accessories like Water Bottle/Food Hopper/Removable floor mesh/Removable Stainless Steel Tray/Data Label Holder etc. Sizes: 450 x 600 x 450 mm.

RABBIT CAGE Stainless Steel :

Stainless Steel Rabbit Breeding/Meeting/Maternity Cage.

Sizes: 36"x18"x18". Complete with all Stainless Steel material.

RABBIT HOLDER :

Box type, Wooden, Glass front, top lid opening, air window, pushing system & neck holding clamp provided. Size 17"x6"x6" with clamp of brass.

RAT METABOLIC CAGE :

Size : 8" x10 ½" x16 ½" (L×B×H). This metabolic cage for **Rats** enables collection of behavioral information, urine and faeces. Made in Stainless Steel 304 quality material & Polycarbonate material. It is also includes feed compartment and water feeding unit.

-

Animal House

Catalog No. 100736

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.