



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Auto Karl Fisher Titrimeter, Digital

Product Image



Description

Auto titration control.
Auto zero burette.
Adjustable timer.
Visual & audible alarm.

Technical Specification:

Display:	Array of 10 LED in line
Sample Size:	1- 50 mg of water
K.F. Dispensing Resolution:	0.05 ml
Indication:	RUN ADD END Display in auto mode
Detection:	Current sensing by electronic control
K.F. Flow:	Highly sensitive solenoid control. Volve for automatic control of K.F. reagent.
Measuring Electrode:	Dual platinum electrode
Titration Cell:	Air tight 200ml glass vessel
End Point Delay:	30 seconds
Buzzer:	Yes
Timer:	Adjustable Timer
Stirrer:	Variable speed magnetic stirrer With On/Off control
Burette:	Pressure filling auto zero burette With 0.05 ml resolution
Power:	230V \pm 10% AC 50Hz
Size:	435 x 305 x 150 mm
Weight:	5 kg (approx.)

Accessories:

Burette 10 ml

Reservoir bottle 500 ml

Rubber bellow for pumping air.

Titration cell 200 ml

Dual platinum electrode.

Teflon coated magnetic stirring capsule.

Nozzle.

Moisture trap with rubber tubing.

Instruction manual.

Particulars**Catalog No.**

Auto Karl Fisher Titrimeter, Digital

JA90110

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.