



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Autoclave, Vertical, Double Walled, Radial Locking

Product Image



Description

Double wall chamber Autoclave. The inner chamber and outer wall made of stainless steel, ring made of stainless steel/ mild steel. The lid is made of stainless steel tightened by radial arms through a steering wheel made of chrome plated mild steel. The instrument is absolutely leak proof. Equipped with pressure gauge, safety valve, double safety valve, steam release valve, rubber gasket, inner container, S.S bucket and mild steel legs for Autoclave. The Autoclave can hold the in between 5psi to 20 psi.

Particulars	Catalog No.
Autoclave, Vertical, Double Walled, Radial Locking, 23 litres (10" x 18")	JA21200
Autoclave, Vertical, Double Walled, Radial Locking 40 litres (12" x 22")	JA21201
Autoclave, Vertical, Double Walled, Radial Locking 52 litres (14" x 22")	JA21202
Autoclave, Vertical, Double Walled, Radial Locking 60 litres (14" x 25")	JA21203
Autoclave, Vertical, Double Walled, Radial Locking 75 litres (16" x 24")	JA21204
Autoclave, Vertical, Double Walled, Radial Locking 95 litres (18" x 24")	JA21205
Autoclave, Vertical, Double Walled, Radial Locking 180 litres (22" x 30")	JA21206

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.