



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Autoclave, Vertical, Stainless Steel, Nut Locking, (Sterilizer Pressure Type)

## Product Image

---



## Description

---

Double wall single chamber Autoclave. The inner chamber, outer wall, lid and ring are made of stainless steel. The instrument is absolutely leak proof. Equipped with pressure gauge, safety valve, double safety valve, steam release valve, vacuum release valve, rubber gasket, inner container, inner tripod stand and mild steel legs for Autoclave. The autoclave can hold the in between 5 psi to 20 psi.

<b>Particulars</b>	<b>Catalog No.</b>
Autoclave, Vertical, Stainless Steel, Nut Locking, 23 litres (10" x 18")	JA21180
Autoclave, Vertical, Stainless Steel, Nut Locking, 40 litres (12" x 22")	JA21181
Autoclave, Vertical, Stainless Steel, Nut Locking, 52 litres (14" x 22")	JA21182
Autoclave, Vertical, Stainless Steel, Nut Locking, 60 litres (14" x 25")	JA21183
Autoclave, Vertical, Stainless Steel, Nut Locking, 75 litres (16" x 24")	JA21184
Autoclave, Vertical, Stainless Steel, Nut Locking, 95 litres (18" x 24")	JA21185
Autoclave, Vertical, Stainless Steel, Nut Locking, 180 litres (22" x 30")	JA21186

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.