



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Automatic Arm Digital Blood Pressure Machine

Product Image



Description

Features :

- Large LCD
- Automatic inflation and deflation

- 3*12 Sets of memory with average
- 3 Minutes Power Off function
- Low battery indicator
- Requires 4x1.5V AA Alkaline batteries
- Cuff size for upper arm circumference between 24-34cm (9.4 to 13.4 inches)
- Date/Time display

Specification:

Measuring methods: Oscillometric

Memory function: 3 x 12 Sets of memory with average digit size: 8.3 x 14.0mm

Display: Digital LCD: 72.8 x 21.8mm

Measuring range: Pressure: 20-280 mmHg Pulse: 40-199 / minute

Sensor accuracy: Pressure: ± 3 mmHg Pulse: $\pm 5\%$

Inflation system: Electro-pneumatic pump

Arm size range: 24 - 34 cm, (9.4-13.4 inches)

Batteries: 1.5V Alkaline (LR6 / AA) x 4

Automatic power-off: Approx. 3 minute after measurement

Reference method for clinical trials: Auscultatory measurement

Weight: Approx. 310g (without batteries)

Storage and transport condition: $-10^{\circ}\text{C} \sim +60^{\circ}\text{C}$, 10% ~ 95% RH

Operating condition: $+10^{\circ}\text{C} \sim +40^{\circ}\text{C}$, 10% ~ 85% RH

Standard cited: EN 1060-1, EN 1060-3, EN 60601-1, EN 60601-1-2, EN 14971 CMDCAS

Catalog No. Particulars

JH10520 Automatic Arm Digital Blood Pressure Machine

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.