



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Automatic Pipette Washer

Product Image



Description

This unit has a large capacity of washing device. Especially recommended for serological, bacteriological or clinical laboratories. It operates on a siphonic principle to send a continuous flow of water over and through the pipettes, at a most normal pressure the

water is completely changed in cycles of about one minute and even at reduced pressure this washer operates efficiently.

It consists of a stainless steel washing cylinder 22" height, 6" diameter mounted on base 10" in diameter and a removable pipette holder which is of 6½" diameter. The cylinder has a ½" inlet with a removable nose nipple for attaching to permanent drain if desired. This unit can accommodate pipettes up to 376mm in length.

Automatic Pipette Washer

Catalog No. 100399

Disclaimer

The Product details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.