



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Automatic Tissue Processor with Thermostat

Product Image



Description

Model with two basket rotor, built in safety device. The equipment is supplied with 2 stainless steel basket rotors, 2 stainless steel baskets, 24 tissue capsules, 10 one litre glass beakers without lip, 2 nos. thermostatically controlled wax bath, 6 nos. of 24 Hr. calibrated metallic timing disc and one cutting plier and dust cover.

Catalog No.	Particulars
JM00530	Automatic Tissue Processor with Thermostat

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.