




Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

BOD Incubator, GMP MODEL

Product Image



Description

Features:

- LCD Display/LED Display.
- Double door design.
- Polished Inner and outer stainless-steel chamber.
- Round angle structure, easy to clean.
- Universal castors for easy movement.
- CFC-free environmental friendly refrigerant.
- Microprocessor PID temperature control: accurate and reliable.
- BOD socket integrated inside for additional small equipment.

Safety function:

- Temperature protection threshold value can be set.
- Auto-off protection and alarm for abnormal situation.
- Parameters can be saved automatically when power off or system error.

Cat No.	Description	Internal Size (WxDxH) mm
101503 A	BOD Incubator, 100L	500 X 350 X 580
101503 B	BOD Incubator, 150L	550 X 401 X 690
101503 C	BOD Incubator, 200L	600 X450 X 750
101503 D	BOD Incubator, 250L	660 X 450 X 850
101503 E	BOD Incubator, 450L	600 X 600 X 1200
101503 F	BOD Incubator, 1000L	800 x 840 x 1500

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.