



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

# BOD Incubator, Stainless Steel

## Product Image



## Description

---

B.O.D. Incubators are suitable for preservation of vaccines, insulin, liver extracts, chemicals etc. and to make biochemical oxygen demand determination.

**Construction :** The low temperature incubator is made of double walled M.S. sheet finished with epoxy based powder coated paint. Inner chamber is made of highly polished stainless steel. Two doors are provided. Inner door is made of transparent acrylic for inspecting specimens outer door is insulated and is fitted with magnetic tape with lock and key. Temperature range from 5° c to 50° c with an accuracy of  $\pm 1^\circ$  c provided with two air circulating fans, high performance compressor, cooling coils, heating elements, digital temperature control and display. It is fitted with a door operated illumination and a caster wheel for easy mobility.

Particulars	Catalog No.
BOD Incubator, Stainless Steel, 825x505x415mm, 6.1 Cubic Feet	JA22200
BOD Incubator, Stainless Steel, 850x600x500mm, 8 Cubic Feet	JA22201
BOD Incubator, Stainless Steel, 880x550x550mm10 Cubic Feet	JA22202
BOD Incubator, Stainless Steel, 900x650x580mm, 12 Cubic Feet	JA22203

**Note:** Custom-made incubators to your specific requirements like size, material, control mechanism, configuration, shelves or special needs available on demand.

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.