



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bar Gel Polyacrylamide Electrophoresis Apparatus Rectangular

Product Image



Description

It is a standard mechanism meant for running 12 gel tubes of 7.5 cm or 14 cm. long. Upper buffer reservoir & lower buffer reservoir is notched from the side to accommodate detachable Cathode & Anode Platinum electrode assembly central compartment helps to create space to accommodate 14 cm. long glass tubes. The elementary unit consist of 3 main detachable compartments, detachable platinum electrode assembly cathode & anode, glass tubes 7.5 cm.(24 Nos.), 14cm.(24 Nos.), silicon rubber grommets(12 Nos.), rubber Bungs(12 Nos.), Polythene Clousers(24 Nos.), Gel casting tube stand, Syringe 5 Ml. with needle, Power cord.

Particulars

Bar Gel Polyacrylamide Electrophoresis Apparatus
Rectangular

Catalog No.

JA20800

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.