



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Bench Top Lab. Centrifuge (Brushless Motor), 5000 r.p.m

Product Image



Description

General purpose Laboratory Centrifuge Machine.
Microprocessor based.
Heavy duty C.R.C. construction finished with powder coating.

Fitted with brushless automatic balanced motor.
 Digital timer 0-99 minutes.
 Digital speed indicator.
 Digital R.C.F. indicator.
 Programmable speed regulator.
 Electronic safety lid lock.
 Imbalance detector.

Special Features

Table top.
 Brushless self balancing motor without carbons.
 High capacity rotors.
 High quality.
 High Performance.

Technical Parameter :

Max. speed: 5000r/min
Max. RCF: 4390xg
Motor: Brushless Motor
Noise: ≤65dB
Power: AC 220V 50Hz 15A
Weight: 35 kg.
Dimensions: 530x420x350 (LxWxH)

Most suitable for medical colleges, clinical labs, where work load is very high.

Details of rotors which can be selected as per customer requirements:

Rotory Type	Rotor Capacity	Speed on Rotor (r.p.m)
Swing Bucket Rotor	4x100ml	5000
Swing Bucket Rotor	4x50ml	5000
Swing Bucket Rotor	8x50ml	4000
Swing Bucket Rotor	32x10ml	4000
Swing Bucket Rotor	32x7ml	4000
Swing Bucket Rotor	48x7ml	4000
Swing Bucket Rotor	48x5ml	4000
Swing Bucket Rotor	72x7ml	4000

Swing Bucket Rotor 4x250ml 4000

Model

Catalog No.

Bench Top Lab. Centrifuge (Brushless Motor), 5000
r.p.m

JA10250

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.