



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Blood Bag Tube Sealer

## Product Image

---



## Description

---

Model : BBS-10

### Construction:

- Portable.
- Body made of CRC sheet duly furnished with powder coating.
- The front of tube Sealer is made of Stainless Steel covering and stainless steel Splash Guard.
- The body is fitted with 4 nibber feet.
- It is also fitted with a Flap for carry the Machine.

### Specifications:

- Microprocessor based control system.
- Very stable and reliable working.
- Working Voltage Range: 220Volts+ 10%. 50 Hz A.C.

### Technical Parameter:

- Sensor: Automatic Senses Blood Bag Tube.
- LED Update : Ready to Seal, sealing in process, sealing completed through LED indication lights.
- Splash Guard: Splash Guard for technician Safety (So Blood don't splash on technician).
- Seal Time: 1.3 Sec. to 1.5 Sec.
- Frequency: RF Fr. 40.68

### Weight and Dimension:

- Weight - 5kg. Approx.
- Dimensions (L•W• H) - 350\*190\*170 mm.

-

Blood Bag Tube Sealer

**Catalog No. 100778**

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.