



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Blood Pressure Machine, Aneroid

Product Image



Description

A compact and reliable Sphygmomanometer supplied with adult velcro fabric cuff and precision manometer with clip for fixing.

Features:

Soft Cotton / Nylon cuff

Emulsoid Bladder & Soft emulsoid Bulb

Leather bag easy to carry
Latex free for safety.
Laser engraved dial face for accuracy
Accurate to: +/- 3mmHg

Catalog No.	Particulars
JH10460	Blood Pressure Machine, Aneroid

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.