



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Blood Pressure Machine, Mercury, Desk Model

Product Image



Description

Mercurial Sphygmomanometer is traditional blood pressure measuring instrument with mercury locking facility. Packed in light weight metal case. Range 0-300 mm Hg.

Features:

Easy to read with clear tube

Accurate reading

Emulsoid Bladder & Soft emulsoid Bulb

Mercury lock facility.

Technical Specification:

Measure Scope: 0-300mmHg

Accuracy: 3mmHg

Glass Tube: 3.5 to 4.0mm

Accessories:

Cotton/Nylon cuff

Latex Bulb and Inflation Valve

Catalog No.	Particulars
--------------------	--------------------

JH10470	Blood Pressure Machine, Mercury, Desk Model
---------	---

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.