



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Blood Pressure Machine, Mercury (Standing Model)

Product Image



Description

Mercurial Sphygmomanometer is traditional blood pressure measuring instrument with mercury locking facility. Packed in lightweight metal case. Range 0-300 mm Hg

Features:

Easy to read with clear tube

Accurate reading.

Emulsoid Bladder & Soft emulsoid Bulb Mercury lock facility. The

Height of pole is adjustable (with range of 105cm to 146cm).

Five wheels with brakes

Technical Specification:

Measure Scope: 0-300mmHg

Accuracy: +/-3mmHg

Glass Tube: 3.5 to 4.0mm

Accessories:

Cotton/Nylon cuff

Latex Bulb and Inflation Valve

Catalog No. Particulars

JH10500 Blood Pressure Machine, Mercury (Standing Model)

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.