



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Blood Pressure Machine, Mercury (Standing Model)

## Product Image

---



## Description

---

Mercurial Sphygmomanometer is traditional blood pressure measuring instrument with mercury locking facility. Packed in lightweight metal case. Range 0-300 mm Hg

**Features:**

Easy to read with clear tube

Accurate reading.

Emulsoid Bladder & Soft emulsoid Bulb Mercury lock facility. The

Height of pole is adjustable (with range of 105cm to 146cm).

Five wheels with brakes

**Technical Specification:**

Measure Scope: 0-300mmHg

Accuracy: +/-3mmHg

Glass Tube: 3.5 to 4.0mm

**Accessories:**

Cotton/Nylon cuff

Latex Bulb and Inflation Valve

**Catalog No.    Particulars**

JH10500	Blood Pressure Machine, Mercury (Standing Model)
---------	--

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.