



**Address:**  
**JAMBU PERSHAD & SONS**  
6275/22 Nicholson Road,  
Ambala Cantt, Haryana,  
INDIA  
Pin: 133001

**Email:**  
sales@japson.com  
japsonambala@yahoo.com

**Website:**  
www.japson.com  
**Phone:**  
+91-171-4006897

---

# Blood Pressure Machine, Mercury Free

## Product Image

---



## Description

---

LCD mercury free sphygmomanometer is a new generation mercury free blood pressure meter which is still made with Krotkoff Sound. It displays results with LED bar graph and completely solves the pollution and the harm of mercury in. LED bar graph sphygmomanometer is an accurate measuring instrument for hospital, clinic and even house using. It does not cause pollution and easy for use, durable. There is a button beside the

LCD screen, if you put the button down, it will be Oscillometric measuring method for the measurement by puffing the bulb manually. It is like digital blood pressure meter, it is used by puff the bulb manually. When the LCD column reaches to about 160-180mmHg, it stops puffing for about 10 seconds and the value displays on to the screen. It is easier than Stethoscopy method.

Catalog No.	Particulars
JH10450	Blood Pressure Machine, Mercury Free

## Disclaimer

---

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.