



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Blood Pressure Machine, Student Model

Product Image



Description

A compact and reliable Sphygmomanometer supplied with adult velcro fabric cuff and precision manometer with clip for fixing.

Features:

Soft Cotton / Nylon cuff

Emulsoid Bladder & Soft emulsoid Bulb

Easy to carry

Latex free for safety

Accurate to: +/- 3mmHg

Please refer to Diagnostic Products Section in Hospital Products Category for detailed specifications and other models.

Catalog No.	Particulars
-------------	-------------

JZ00290	Blood Pressure Machine, Student Model
---------	---------------------------------------

A quality product from JAPSON ideal for use in School Laboratories in developing countries of Asia, Africa and South America.

Standard specifications for school laboratory use. Ideal mix of Quality and Economy.

Other configurations and superior quality available for specific use or for developed countries.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.