



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Blood Roller Mixer

Product Image



Description

Catalog No : 100089

Blood Roller Mixer - Model 43

Four Roller design with small footprint

SALIENT FEATURES

- ✓ Rocking and rolling action for complete Mixing
- ✓ Four roller design with small footprint
- ✓ Designed for continuous quiet operation
- ✓ Ambient temperature range +4°C to +60°C

TECHNICAL SPECIFICATIONS

✓ ROLLER SIZE	230mm
✓ CONTROLS	Analog
✓ NO. OF ROLLERS	4
✓ SPEED	33rpm (Approx.)
✓ MAXIMUM LOAD	5 test tubes
✓ TIMER	No
✓ POWER SUPPLY	230V \pm 10% AC 50Hz.

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.