



Address:
JAMBU PERSHAD & SONS
6275/22 Nicholson Road,
Ambala Cantt, Haryana,
INDIA
Pin: 133001

Email:
sales@japson.com
japsonambala@yahoo.com

Website:
www.japson.com
Phone:
+91-171-4006897

Box for Micro Centrifuge Tubes

Product Image



Description

Material: Polypropylene

The tight fit cover ensures samples won't fall if box is reversed. This box is autoclavable at 121°C & freezable at -80°C.

Catalog No.

JW08641

JW08642

Capacity (ml)

0.5 ml/75 places

1.5 ml/75 places

JW08643

2 ml/75 places

Disclaimer

The Products details given on this page are indicative in nature and JAPSON reserves the right to change them without prior notice. Buyer is also requested to re-check the specifications and other features of product at the time of order as product development is a continuous process and minor modifications may be made to design based on latest availability, process and design.